

NSK's Observation on the research paper

"In vitro study of anti-suck-back ability by themselves on new high-speed air turbine handpieces"

Written by Toshiko OZAWA, Masako NAKANO, and Takashi ARAI
 Published by Dental Materials Journal (2010)



- **The Test Method & Result of The Paper**
- **Our Observation on The Paper**
- **Clean Head System** (NSK Anti-Suckback system)
- **Our Experiment**
- **Result of Our Experiment**




Fig. 1
Complete test setup : Measuring suckback pressure of the whole handpiece brand.




Fig. 2
Face test setup : Measuring suckback pressure of each face of the handpiece.



Fig. 3
Test setup to verify suckback Mist condition

※ These figures are extracted from the paper

One of the test result of the paper

Maximum water height Mean (SD)mm, n=5

Handpiece model	Whole head	Front face	Back face
A	426.0 (5.5)	358.0 (4.5)	219.0 (1.2)
B	31.4 (0.5) *	34.6 (0.5)	<0

* no significant difference (p>0.01)

According to this paper, Sample ID A, or NSK Ti-Max X700L, showed the highest water height in the table "Maximum negative pressure of turbine handpieces measured as the maximum water height." Table 2, p652.



Fig. 1
Complete test setup : Measuring suckback pressure of the whole handpiece brand.




Fig. 2
Face test setup : Measuring suckback pressure of each face of the handpiece.





Fig. 3
Test setup to verify suckback Mist condition

All three test method were done in sealed condition.

During normal clinical treatment, there are no circumstances where turbine is operated under such sealed state.


NSK's Anti Suckback system called "Clean Head System" uses the pressure difference between inside and outside of the turbine to block extraneous substance going into the head of handpiece. This test method of the paper will create misleading perception to the effectiveness of the Anti-Suckback function of our Clean Head System.

During Rotation



During the rotation, contamination object will not enter to the inside of head. The inside air pressure is higher than outside, therefore the air flows from inside to outside.

Right After Stopping



After the air supply stopped, cartridge still rotate with momentum. "parts A" will rotate together with rotor and the contamination object in mouth will be flown out from head by centrifugal force.

Objective To measure the effectiveness of NSK Clean Head System

Testing Procedure

Submerge tip of head 1mm below the red Iodine Solution



1. Operate Turbine for 10 sec



2. Stop Turbine for 10 sec



Repeat
10 times



After experiment, red iodine traces can be found on the following:




	NSK		K		S	M		W		B
	Z900L	M900L	K1	K2	S	M1	M2	W1	W1	B
A	○	○	○	×	×	×	×	×	×	×
B	○	○	○	×	×	×	×	×	×	×

× = Red iodine traces found

NSK Z900L & M900L and K1 did not suck back red iodine at all. However, the other companies' turbines sucked back not only to the rotor, but also to the hose.

Suckback Test Results: weight of suckback water



The amounts of residue sucked back by NSK Z900L & M900L and K1 were very minimal.
However, the ones of the other companies' turbines were substantially larger.



CREATE IT.